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

www.reh4mat.com

SPLINT 1X USER MANUAL - 1X

Orthopedic splint 1X is a advanced product that allows you to adjust the range of motion of your joint – flexion, extension or immobilization in the required position. The splint is protected against accidental or unauthorized change of the joint angle.

ATTENTION:

The change of ROM should be only performed by the doctor or physiotherapist.

1. FLEXION:

The flexion of the knee joint is adjusted with the lower slider (1). To change the flexion angle, press the round button in the central part of the slider and move it to the appropriate flexion value. After setting the button in a right place, release it. The button will pop out and lock the slider against the accidental movement by the user.

The splint is equipped with the following angles of joint flexion:

full mobility, 90°, 70°, 40°, 30°.

-10 > 30°

2. EXTENSION:

The extension angle is adjusted using the upper slider (2). To change the extension angle, press the round button in the central part of the slider and move it to the appropriate extension value. After setting the button in a right place, release it. The button will pop out and lock the slider against the accidental movement by the user.

The splint is equipped with the following angles of joint extension:

full mobility, 30°, 10°, 0°, -10°.

>10°

3. CONSTANT JOINT LOCK:

It is possible to permanently lock in the required angle of the joint, what is very useful when it is necessary to immobilize a limb. The angle is locked in a given setting using a slider (3) in the upper part of the splint. To set the angle of the joint lock, press the round button in the central part of the slider and, holding it, flex the splint in the axis of its rotation. When the appropriate angle value is opposite the slider, move it clockwise and release the round button. The slider will lock automatically when the round button will pop up and the metal part slides into the appropriate value.

The splint can be locked in the following angles:
full mobility, 90°, 80°, 70°, 60°, 50°, 40°, 30°, 20°, 10°, 0°, -10°.

The skin is protected with special soft cushions, so using the splint is comfortable.

>40°<

ATTENTION:

If the brace is equipped with two side splints, the adjustment should always be made separately on each of the splints, setting the range of motion at the same angle of flexion, extension or lock.

4. SPLINT LENGTH ADJUSTMENT:

Each 1X splint is equipped with two upper and two lower plastic pads, through which the circumferential straps run. The distance between pads can be adjusted by pressing the oval button (4), moving it to the appropriate position and releasing the button. When the oval button pops up, the pad locks in given place.


A. Symmetrical setting

B. Asymmetrical setting

ATTENTION:

The pads on the upper part of the splint can be set in different position than the pads on the lower arm. The upper part of the splint may have a different length than the lower part - however, the rule is that if there are two side splints in a given orthosis, the upper arms in both splints must be identical and the lower arms must be identical.

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