


Orteza została wyposażona w szynę boczną, która posiada **blokadę kąta zgięcia (A)** w zakresie: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120° oraz **blokadę wyprostu (B)** w zakresie: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°.




1. Ustaw ortezę w przybliżonym położeniu, w którym chcesz ograniczyć **zgięcie/wyprost**
2. Wkręć wkręt/wkręty w otwór oznaczony odpowiednią wartością kąta i sprawdź poprawność działania oraz zakres ruchów zegara (**rys.1**)
3. **(A) - otwory oznaczone na powyższym zdjęciu obszarem (A) odpowiadają za zgięcie (rys.2)**  
**(B) - otwory oznaczone na powyższym zdjęciu obszarem (B) odpowiadają za wyprost (rys.2)**

**ZAKRES RUCHÓW I JEGO WARTOŚCI KĄTOWE USTALA LEKARZ PROWADZĄCY!**

**UWAGA:** Jeżeli orteza wyposażona jest w dwie szyny boczne, regulacji należy dokonać osobno na obydwu szynach, ustawiając ten sam kąt zgięcia lub wyprostu. Szyny boczne nie wymagają konserwacji

## HOW TO ADJUST THE FLEXION AND EXTENSION ANGLE IN 1R SPLINT - INSTRUCTION

1R SPLINT is equipped with the following ROM adjustment **flexion (A)** in the scope of: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120° and **the extension (B)** in the scope of: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°.



1. Position the brace in the approximate position where you want to limit the **flexion/extension**
2. Drive the screw/bolts into the hole marked with the appropriate angle value and check the correct operation and the range of clock movements (**fig.1**)
3. **(A) - The holes indicated in the picture as (A) are responsible for flexion (fig.2)**  
**(B) - The holes indicated in the picture as (B) are responsible for extension (fig.2)**

**THE ANGLE OF FLEXION AND EXTENSION IS SPECIFIED BY THE DOCTOR OR PHYSIOTHERAPIST!**

**ATTENTION:** If the brace is equipped with 2 side splints, adjustment is done on each side splint separately setting the same angle of flexion and extension. Side splints do not require maintenance



REH4MAT

Ul. Piaski 47, 36-060 Głogów Małopolski | Polska

tel. +48 16 62 142 20

Fax +48 16 62 21 42 13

e-mail: [biuro@reh4mat.com](mailto:biuro@reh4mat.com)

[www.probraces.com](http://www.probraces.com)

[www.reh4mat.com](http://www.reh4mat.com)